

Baby Infant Warmer

AMB-800


Features

- ✓ Microprocessor based servo controlled temperature system
- ✓ Three modes: pre-warm, manual or automatic
- ✓ Audible and visual Alarm functions :
- ✓ Power failure, temperature deviation, temperature sensor failure, skin temperature failure, over temperature
- ✓ 8"LCD colorful touch screen
- ✓ Carbon steel heater more effective
- ✓ Special design inclined heater head is more convenience for X ray test image
- ✓ Digital skin temperature sensor,the data is accurate and stable
- ✓ Alloyed bracket,fastness
- ✓ Integral timer to record therapy time
- ✓ Bigger mattress
- ✓ Skin temperature trend for 2 hours
- ✓ Electric infant bed tilt angle is $\pm 12^\circ$
- ✓ Four transparent protector folded with 2 infusion seals
- ✓ Convenient X-ray cassette tray under infant bed
- ✓ Electric height adjusted suit for operation
- ✓ With observe light,oxygen supply system and LED phototherapy
- ✓ With RS-232 connector,2 trays and 2 drawers



TECHNICAL DATA

Power supply	: AC110/220V 60/50Hz
Power input	: 1100VA
Operating conditions	Environment temperature : 18°C ~ 30°C
	Environment relative humidity : 30 ~ 75
	Atmospheric Pressure : 700 ~ 1060hpa
	Environment air velocity of flow : 1/2 P V
Skin Temperature control range	: 32°C ~ 38°C
Temperature control accuracy	: ≤0.5°C
Skin temperature sensor accuracy	: ±0.3°C
Mattress temperature uniformity	: 1/2 •
Infant Bed tilt angle	: ±12°
Warm-up time(from 25°C)	: 1/2 P IQ
LED observation light	: 3W
Distance from heater to mattress	: 80cm
Mattress size	: 78cm X 64cm
Timer range	: 0 ~ 99h59min
Height adjustable range	: 0 ~ 20cm
LED phototherapy unit 36cm	: 620-2800μw/c
LED phototherapy unit 50cm	: 400-1620μw/c
Oxygen cylinder	: 6L X 2
Oxygen flow meter	: 0-10L/min

**Infant
Resuscitator**

**Low pressure
aspirator**

**Air&Oxygen
Blender**

STANDARD CONFIGURATION

Main body (including the Radiant source, Control system, Infant bed, Bracket), X-ray cassette, Observe light, I.V. pole, Skin temperature sensor, Tray, Mattress, Transparent protector, Castors, Drawers, Oxygen supply system, LED phototherapy unit, Electric height

Gel mattress

PACKAGE

Each unit is packed one plywood case
Size: 118*90*87cm

N.W. : 100Kg
G.W. : 130Kg

TRANSPORT & STORAGE

Environment temperature -40°C ~ +55°C
Environment relative humidity ≤95
Atmospheric Pressure 500 ~ 1060hpa



55 The Avenue
London
E30 5AD
www.astramed.co.uk

