

Compression Therapy

CT-01

Features

- ✓ Skip care and intensive care optional to optimize treatment effect
- ✓ Adjustable intensity level and free time setting from 10, 20, 30mins
- ✓ Treatment process controlled by computer program automatically and accurately
- ✓ Color LED touch screen for more convenient use
- ✓ Indicator lights will be on while related air chamber working
- ✓ Handle pick-up design for more convenient carry
- ✓ Could connect with different sleeves or cuffs separately or simultaneously




**Overlapping to
prevent blood
dot backflow**

**Normal
Sleeve for economic
choosing**



3 Massage Modes



Specification

Machine Size:	210*256*227mm
Machine Weight:	about 1.5KGS
Pressure Range:	20-200mmHg
Cuffs Color:	Black, Grey, or Customized Color
Input Voltage:	100-240V 50Hz/60Hz
Time Setting:	10, 20, 30Mins
Pressure Level:	10
Power Consumption:	18w

Legs Cuffs Dimension

SIZE	A	B	C	D	E	F	G
M (non-overlapping) (91*65cm)	65cm	51cm	41cm	33cm	77cm	36cm	27cm
M (non-overlapping) (91*65cm)	74cm	55cm	41cm	34cm	80cm	39cm	30cm
M (non-overlapping) (91*65cm)	70cm	60cm	54cm	46cm	86cm	43cm	35cm
M (non-overlapping) (91*65cm)	75cm	65cm	54cm	46cm	100cm	53cm	35cm

Extension Zipper of width 10cm to widen legs cuffs for Obesity Use (optional Accessory)